



Version: 1.0
Initial preparation date: 2023.02.23
Web:www.yacoo.com.cn
E-mail:sales@yacoo.com.cn

Material Safety Data Sheet

SECTION 1:CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: TAPS

Company: Suzhou YACOO Science Co., Ltd.

Address: No.128,FangZhou Road,Suzhou Industrial Park,China

Tel: 0512-87182055

Fax: 0512-87182056

SECTION 2: Hazards identification

Summary of emergency

crystalline white After inhalation: fresh air. In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/ shower. After eye contact: rinse out with plenty of water., Remove contact lenses. After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell. Vapors are heavier than air and may spread along floors. Forms explosive mixtures with air on intense heating.

2.1GHS Classification

Not a hazardous substance or mixture according to the Globally Harmonized System (GHS).

2.2GHS Label elements, including precautionary statements

Not a hazardous substance or mixture according to the Globally Harmonized System (GHS).

2.3Physical and chemical hazards

Referring to current information, no physical or chemical hazard.

2.4Health hazards

Referring to current information, no health hazard.

2.5Environmental hazards

Referring to current information, no environmental hazard.

2.6Other hazards - none

SECTION 3: Composition/information on ingredients

Substance / Mixture: Substance

3.1Substances

Synonyms: [(2-Hydroxy-1,1-bis(hydroxymethyl)ethyl)amino]-1-propanesulfonic acid

N-[Tris(hydroxymethyl)methyl]-3-aminopropanesulfonic acid

Formula: C7H17NO6S



Molecular weight: 243.28 g/mol

CAS-No.: 29915-38-6

EC-No.: 249-954-1

No components need to be disclosed according to the applicable regulations.

SECTION 4: First aid measures

4.1 Description of first-aid measures If inhaled

After inhalation: fresh air.

In case of skin contact

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

In case of eye contact

After eye contact: rinse out with plenty of water. Remove contact lenses.

If swallowed

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of any immediate medical attention and special treatment needed

No data available

4.4 Notes to physician

No data available

SECTION 5: Firefighting measures

5.1 Extinguishing media Unsuitable extinguishing media

For this substance/mixture no limitations of extinguishing agents are given.

5.2 Special hazards arising from the substance or mixture

Carbon oxides Nitrogen oxides (NOx) Sulfur oxides

Vapors are heavier than air and may spread along floors. Forms explosive mixtures with air on intense heating.

5.3 Advice for firefighters

In the event of fire, wear self-contained breathing apparatus.

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures Advice for non-emergency personnel: Avoid inhalation of dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

For personal protection see section 8.

6.2 Environmental precautions

Do not let product enter drains.

6.3 Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills. Observe possible material restrictions (see sections 7 and 10). Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

6.4 Reference to other sections
For disposal see section 13.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities Storage conditions

Tightly closed. Dry.

Storage class

Storage class (TRGS 510): 13: Non Combustible Solids

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Ingredients with workplace control parameters

We are not aware of any national exposure limit.

8.2 Exposure controls

Appropriate engineering controls

Change contaminated clothing. Wash hands after working with substance.

Personal protective equipment Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU). Safety glasses

Skin protection

This recommendation applies only to the product stated in the safety data sheet, supplied by us and for the designated use. When dissolving in or mixing with other substances and under conditions deviating from those stated in EN374 please contact the supplier of CE-approved gloves (e.g. KCL GmbH, D-36124 Eichenzell, Internet: www.kcl.de).

Full contact

Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

Material tested: KCL 741 Dermatril® L

Splash contact Material: Nitrile rubber

Minimum layer thickness: 0.11 mm

Break through time: 480 min

Material tested: KCL 741 Dermatril® L

Respiratory protection

required when dusts are generated.

Our recommendations on filtering respiratory protection are based on the following standards: DIN EN 143, DIN 14387 and other accompanying standards relating to the used respiratory protection system.

Control of environmental exposure

Do not let product enter drains.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

- a) Physical state crystalline
 - b) Color white
 - c) Odor No data available
 - d) Melting point/freezing point Melting point/range: 230 - 235 ° C - dec.
 - e) Initial boiling point and boiling range No data available
 - f) Flammability (solid, gas) No data available
 - g) Upper/lower flammability or explosive limits No data available
 - h) Flash point > 110 ° C - closed cup
 - i) Autoignition temperature No data available
 - j) Decomposition temperature No data available
 - k) pH No data available
 - l) Viscosity
Viscosity, kinematic: No data available
Viscosity, dynamic: No data available
 - m) Water solubility ca. 121.6 g/l at 20 ° C
 - n) Partition coefficient: n-octanol/water No data available
 - o) Vapor pressure No data available
 - p) Density No data available
Relative density No data available
 - q) Relative vapor density No data available
 - r) Particle characteristics No data available
 - s) Explosive properties No data available
 - t) Oxidizing properties No data available
- ## 9.2 Other safety information
- No data available

SECTION 10: Stability and reactivity

10.1 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

10.2 Possibility of hazardous reactions

No data available

10.3 Conditions to avoid

Strong heating.

10.4 Incompatible materials

Strong oxidizing agents

10.5 Hazardous decomposition products

In the event of fire: see section 5

SECTION 11: Toxicological information

11.1 Information on toxicological effects Acute toxicity

Oral: No data available Inhalation: No data available Dermal: No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitization

No data available
Germ cell mutagenicity
No data available
Carcinogenicity
No data available
Reproductive toxicity
No data available
Specific target organ toxicity - single exposure
No data available
Specific target organ toxicity - repeated exposure
No data available
Aspiration hazard
No data available
11.2 Additional Information
To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

SECTION 12: Ecological information

12.1 Toxicity
No data available
12.2 Persistence and degradability
No data available
12.3 Bioaccumulative potential
No data available
12.4 Mobility in soil
No data available
12.5 Results of PBT and vPvB assessment
PBT/vPvB assessment not available as chemical safety assessment not required/not conducted
12.6 Endocrine disrupting properties
No data available

12.7 Other adverse effects
No data available

SECTION 13: Disposal considerations

13.1 Waste treatment methods Product
Offer surplus and non-recyclable solutions to a licensed disposal company.

SECTION 14: Transport information

14.1 UN number
ADR/RID: -IMDG: -IATA-DGR: -
14.2 UN proper shipping name ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA-DGR: Not dangerous goods
14.3 Transport hazard class(es)

ADR/RID: -IMDG: -IATA-DGR: -

14.4Packaging group

ADR/RID: -IMDG: -IATA-DGR: -

14.5Environmental hazards

ADR/RID: noIMDG Marine pollutant: noIATA-DGR: no

14.6Special precautions for user

14.7Incompatible materials

Strong oxidizing agents

Further information

Not classified as dangerous in the meaning of transport regulations.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulatory information

Measures on the Environmental Administration of New Chemical Substances Registration

Registration/Notification number: B1A222211757

Downstream users need to comply with the conditions of safe use of the chemical, understand the environmental and health hazard and risk management measures identified on the SDS as well as the local/national regulations concerning the chemical.

Other regulations

Please pay attention on the waste treatment should also comply with local regulations requirement.

SECTION 16:OTHER INFORMATION

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no way shall the company be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if the company has been advised of the possibility of such damages.